

MK-2 Pharm E-Cell® Stack

E-Cell®

E-Cell®, the industry standard for chemical-free EDI technology. E-Cell®'s robust, modular electrodeionization (EDI) technology is cost-effective for all flow rates.

No regeneration chemicals are needed, no hazardous waste stream is produced, operation is simple and continuous, and facility requirements are reduced.

MK-2 Pharm



E-Cell®

MK-2 Pharm

Process Design Parameters		
Parameter	US Units	Metric Units
System Product Flow	7 to 18 gpm	1.59 to 4.09 m ³ /h
Nominal Recovery	80 to 95%	80 to 95%
Temperature	40 to 100° F	4.4 to 38° C
Inlet Pressure	50 to 100 psig	3.4 to 6.9 bar
Feed to Product delta P	25 to 43 psi	1.7 - 3.0 bar
Dimensions	12" W x 19" D x 24" H	30cm W x 48cm D x 61cm H
Product Pipe Material	PP	PP

Operating Conditions			
Condition		US Units	Metric Units
Electrical	Maximum	4.5 amps / Stack @ 600V DC	4.5 amps / Stack @ 600V DC
Product Outlet	Flow Range	7 to 18 gpm	1.59 to 4.09 m ³ /h
	Quality	>10 MOhm.cm	>10 MOhm.cm
	Pressure Drop	25 to 43 psi	1.7 to 3.0 bar
	Temperature Rise	4.1° F maximum	2.4° C maximum
Electrolyte Outlet	Flow	0.25 to 0.40 gpm / Stack to Drain	0.95 to 1.51 lpm/Stack to Drain
	pH	7.0 to 9.0	7.0 to 9.0
Concentrate Bleed	Flow	Determined by Recovery Rate	Determined by Recovery Rate
Concentrate + Electrolyte Inlet	Maximum Flow	5.4 gpm per Stack	1.23 m ³ /h per Stack
	Pressure	10 psi < Feed Water Pressure	0.7 bar < Feed Water Pressure
	Concentrate Conductivity	50 to 800 uS/cm	50 to 800 uS/cm
Concentrate	Flow	Determined by Recovery Rate	Determined by Recovery Rate
Make-up	Water Quality	Same as Feed Inlet	Same as Feed Inlet

Note: The feed water to the E-CELL® system must be RO permeate or equivalent.



E-Cell®

E-Cell Corporation

52 Royal Road, Guelph, Ontario, Canada N1H 1G3
 Tel: +1 (519) 836-2260 Fax: +1 (519) 836-0982
 Toll Free: 1 (888) 386-5132 (USA and Canada)
 Website: www.ecell.com E-mail: ecell@ecell.com